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## Instructions for Use

### CelliMax<sup>®</sup> MetaCHO 培养基

### CelliMax<sup>®</sup> MetaCHO Medium

Cat. No.: CPDP050 (Powder)  
CPDP051 (Liquid)

Issued on January 1, 2026

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##### Introduction

- 适用于 CHO 细胞的培养  
For CHO cell culture
- 化学成分确定  
CD (Chemically Defined)
- 无动物源性成分  
ADCF (Animal Derived Component Free)
- 不含血清  
Serum free
- 不含有酚红  
Without Phenol Red
- 不含有 L-谷氨酰胺  
Without L-Glutamine
- 不含有水解物和重组蛋白  
Without Hydrolysates and Recombinant Protein
- 含有 5.6g/L 葡萄糖  
With 5.6g/L Glucose

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##### Storage

1. 干粉培养基密封保存在 2~8°C 的避光条件下。  
Powder medium should be sealed at 2 to 8°C, away from light.
2. 液体培养基密封保存在 2~8°C 的避光条件下。  
Liquid medium should be sealed dry at 2 to 8°C, away from light.

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##### Shelf Life\*

1. 干粉培养基未开封情况下有效期 2 年, 开封后请尽快使用。详细信息请见标签。  
Powder medium is valid for 2 years without opening. Please use as soon as possible after opening. Please refer to the label for details.
2. 液体培养基未开封情况下有效期 1 年, 开封后请尽快使用。详细信息请见标签。  
Liquid medium is valid for 1 year without opening. Please use as soon as possible after opening. Please refer to the label for details.

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\*The shelf life is based on the general testing method of CelluPro. It is suggested to study on validity for different uses or projects.

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##### Preparation Recommendations

- ◆ 洁净的配液容器。  
Clean vessel.
- ◆ 高纯度水, 如注射用水、超纯水等。  
High purity water (HPW) such as water for injection (WFI) and ultrapure water (UPW).
- ◆ 搅拌效率充足的配液系统。  
Liquid mixing system with sufficient mixing efficiency.
- ◆ 校准的 pH 计 (建议校准点: 4.01、7.00、

9.21)。

Calibrated pH meter (Recommended calibration points: 4.01, 7.00, 9.21).

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Calibrated osmometer (Recommended calibration points: 0, 300, and 700mOsmol/kg).

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Small-scale preparation less than 3 liters is not recommended to prevent errors.

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All utensils or consumables in contact with the medium shall be mycoplasma and endotoxin contamination risk free for product quality.

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Filtration system and liquid storage system in contact with the filtered liquid medium should be kept sterile.

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Pre filtration at 0.8 $\mu$ m or 0.45 $\mu$ m is recommended to improve filtration speed and filtration load.

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Dissolution time varies with preparation scales and mixing efficiency. All the stirring time herein is a reference for small-scale preparation. For larger preparation scale or lower mixing efficiency, each mixing time needs to extend appropriately to ensure sufficient mixing and dissolution effect.

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Please shorten the total medium preparation time as much as possible under ideal mixing effect and powder dissolution effect, for long preparation time will pose a higher risk of bioburden and endotoxin contamination.

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**Q.S. to final volume:** q.s. the preparation solution to the final volume ( $V_{\text{final}}$ ), which is recommended for small-scale preparation.

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**Q.S. to final weight:** q.s. the preparation solution to the final weight ( $m_{\text{final}} = \rho \times V_{\text{theoretical}}$ ), which is recommended for large-scale preparation.

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Recommended density of CelliMax<sup>®</sup> MetaCHO Medium:  **$\rho = 1.008 \text{ kg/L}$** . (Applied to q.s. to the final weight\*)

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\*Take the method of q.s. to the final weight to prepare 1000L ( $V_{\text{theoretical}}$ ) scale of CelliMax<sup>®</sup> MetaCHO Medium as an example: the final weight ( $m_{\text{final}}$ ) is calculated as:  $m_{\text{final}} = \rho \times V_{\text{theoretical}} = 1.008 \text{ kg/L} \times 1000 \text{ L} = 1008 \text{ kg}$ . Then q.s. the medium solution to the final weight of 1008kg, at this time, the final volume of the preparation solution is 1000L.

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Before the preparation of the medium, please

determine the most suitable q.s. method according to the target preparation volume and the available q.s. device (standard volume measuring instruments or weighing equipment).

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CelliMax® MetaCHO Medium is a protein-free CD basal medium. If the original medium with protein components (e.g. insulin and IGF) is replaced with CelliMax® MetaCHO Medium, corresponding substances can be added to shorten the cell adaptation time.

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\* Except for process changes, other changes will not be notified.

©2026 CelluPro, Ltd. All rights reserved including data, graphics and images. CelliMax, CelluPro and Helix logo are registered trademarks of CelluPro, Ltd.