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Instructions for Use

CelliMax[®] MultiCHO 培养基

CelliMax[®] MultiCHO Medium

Cat. No.: CPDP077 (Powder)
CPDP095 (Liquid)

Issued on January 1, 2026

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Introduction

- 适用于 CHO 细胞的培养
For CHO cell culture
- 化学成分确定
CD (Chemically Defined)
- 无动物源性成分
ADCF (Animal Derived Component Free)
- 不含血清
Serum free
- 不含有酚红
Without Phenol Red
- 不含有 L-谷氨酰胺
Without L-Glutamine
- 不含有水解物和重组蛋白
Without Hydrolysates and Recombinant Protein
- 含有 5.6g/L 葡萄糖
With 5.6g/L Glucose

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Storage

1. 干粉培养基密封保存在 2~8°C 的避光条件下。
Powder medium should be sealed at 2 to 8°C, away from light.
2. 液体培养基密封保存在 2~8°C 的避光条件下。
Liquid medium should be sealed at 2 to 8°C, away from light.

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Shelf Life*

1. 干粉培养基未开封情况下有效期 2 年, 开封后请尽快使用。详细信息请见标签。
Powder medium is valid for 2 years without opening. Please use as soon as possible after opening. Please refer to the label for details.
2. 液体培养基未开封情况下有效期 1 年, 开封后请尽快使用。详细信息请见标签。
Liquid medium is valid for 1 year without opening. Please use as soon as possible after opening. Please refer to the label for details.

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*The shelf life is based on the general testing method of CelluPro. It is suggested to study on validity for different uses or projects.

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Preparation Recommendations

- ◆ 洁净的配液容器。
Clean vessel.
- ◆ 高纯度水, 如注射用水、超纯水等。
High purity water (HPW) such as water for injection (WFI) and ultrapure water (UPW).
- ◆ 搅拌效率充足的配液系统。
Liquid mixing system with sufficient mixing efficiency.
- ◆ 校准的 pH 计 (建议校准点: 4.01、7.00、9.21)。

Calibrated pH meter (Recommended calibration points: 4.01, 7.00, 9.21).

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Please shorten the total medium preparation time as much as possible under ideal mixing effect and powder dissolution effect, for long preparation time will pose a higher risk of bioburden and endotoxin contamination.

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Q.S. to final volume: q.s. the preparation solution to the final volume (V_{final}), which is recommended for small-scale preparation.

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Q.S. to final weight: q.s. the preparation solution to the final weight ($m_{\text{final}} = \rho \times V_{\text{theoretical}}$), which is recommended for large-scale preparation.

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Recommended density of CelliMax[®] MultiCHO Medium: **$\rho = 1.008 \text{ kg/L}$** . (Applied to q.s. to the final weight*)

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*Take the method of q.s. to the final weight to prepare 1000L ($V_{\text{theoretical}}$) scale of CelliMax[®] MultiCHO Medium as an example: the final weight (m_{final}) is calculated as: $m_{\text{final}} = \rho \times V_{\text{theoretical}} = 1.008 \text{ kg/L} \times 1000 \text{ L} = 1008 \text{ kg}$. Then q.s. the medium solution to the final weight of 1008kg, at this time, the final volume of the preparation solution is 1000L.

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Before the preparation of the medium, please determine the most suitable q.s. method according to the target preparation volume and the available q.s. device (standard volume measuring instruments or weighing equipment).

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**The pH will gradually increase by the mixing time after dissolution of sodium bicarbonate powder, and the mixing time should be controlled reasonably.*

6. 搅拌时间结束后, 检测培养基溶液 pH。**缓慢**加入 6mol/L 氢氧化钠溶液, 调整培养基溶液 pH 至 **7.00~7.30***。
After stirring, test the pH of the medium solution, and **slowly** add 6mol/L NaOH solution to adjust the pH to **7.00 to 7.30***.
**pH 调整时应注意控制氢氧化钠溶液的使用量, 避免过度调整导致培养基最终渗透压超出标准范围。
The addition amount of NaOH solution shall be controlled to adjust the pH, avoiding out-of-specification osmolality caused by excessive adjustment.
7. 根据配制开始前确定的定容方式, 使用高纯度水定容至终体积或终重量*, 继续搅拌 **10~30 分钟**。
Add HPW to the final volume or final weight* according to the determined q.s. method before preparation and stir for **10 to 30 minutes**.
**培养基定容的具体操作和注意事项详见本说明书中“定容方式”部分的相关描述。
Specific q.s. operations and precautions are described in the “Q.S. Method” section.
8. 此时检测渗透压应为 **290~330mOsmol/kg**, 浊度应 **<3.0NTU**。
At this time, expected osmolality is **290 to 330mOsmol/kg**, and expected turbidity is **< 3.0NTU**.
9. 使用 **0.22µm 或 0.2µm** 除菌滤器*过滤除菌。
Sterilize with **0.22µm or 0.2µm** sterile filter*.
**建议使用低结合类型滤膜, 如聚偏二氟乙烯 PVDF、聚醚砜 PES 或醋酸纤维素 CA。
It is recommended to use a low-binding filter membrane type, such as PVDF, PES or CA.
10. 立即使用, 或密封后在 **2 至 8°C** 条件下**避光**保存。
Use the prepared medium immediately or seal and store at **2 to 8°C away from light**.

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*Please visit www.cellupro.com for more information or contact sales@cellupro.cn.

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* Except for process changes, other changes will not be notified.

©2026 CelluPro, Ltd. All rights reserved including data, graphics and images.
CelliMax, CelluPro and Helix logo are registered trademarks of CelluPro, Ltd.