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## Instructions for Use

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#### Storage

1. 干粉培养基密封保存在 2~8°C 的避光条件下。  
Powder medium should be sealed at 2 to 8°C, away from light.
2. 液体培养基密封保存在 2~8°C 的避光条件下。  
Liquid medium should be sealed at 2 to 8°C, away from light.

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#### Shelf Life\*

1. 干粉培养基未开封情况下有效期 2 年, 开封后请尽快使用。详细信息请见标签。  
Powder medium is valid for 2 years without opening. Please use as soon as possible after opening. Please refer to the label for details.
2. 液体培养基未开封情况下有效期 6 个月, 开封后请尽快使用。详细信息请见标签。  
Liquid medium is valid for 6 months without opening. Please use as soon as possible after opening. Please refer to the label for details.

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\*The shelf life is based on the general testing method of CelluPro. It is suggested to study on validity for different uses or projects.

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#### Preparation Recommendations

- ◆ 洁净的配液容器。  
Clean vessel.
- ◆ 高纯度水, 如注射用水、超纯水等。  
High purity water (HPW) such as water for injection (WFI) and ultrapure water (UPW).
- ◆ 搅拌效率充足的配液系统。  
Liquid mixing system with sufficient mixing efficiency.

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- ◆ 过长的配制时间会带来更高的微生物负荷及细菌内毒素污染的风险，因此在搅拌系统效率和干粉溶解效果理想的条件下，请尽可能缩短培养基配制总时长。  
Please shorten the total medium preparation time as much as possible under ideal mixing effect and powder dissolution effect, for long preparation time will pose a higher risk of bioburden and endotoxin contamination.

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the target preparation volume and the available q.s. device (standard volume measuring instruments or weighing equipment).

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*\*The preparation temperature is for reference only. Very low temperature will reduce dissolution efficiency. The water temperature is recommended to be adjusted to 25°C during medium preparation to enhance the solubility of the materials, for the perfusion medium is nutrient-rich.*

2. 开始搅拌。搅拌速度调整为可使粉末快速地完全浸没在溶液中, 但**不会产生大量气泡**。  
Start stirring. Mix at a suitable speed until completely dissolved **without many bubbles generated**.
3. **缓慢**加入培养基干粉 **34.52g/L\***, 避免大块干粉直接加入水中, 搅拌**不少于 30 分钟**, 此时溶液**仍有浑浊**为正常现象。  
While stirring, **slowly add 34.52g/L\*** powder medium to the vessel, avoiding formation of clumps, and mix **at least 30 minutes**. The solution will **remain cloudy** in this step.  
*\*本文中涉及的所有物料重量, 应尽可能保证最小的误差(1%以内)。若误差无法避免, 应结合实际情况进行研究。*  
*\*The weight error should be kept in minimum (within 1%) for all materials involved herein. If the error is unavoidable, it should be studied based on the actual situation.*
4. **缓慢**加入 6mol/L 氢氧化钠溶液\*, 调整培养基溶液 pH 至 **6.60~6.80\***, 继续搅拌至少 15 分钟, 此时溶液应逐渐澄清。  
**Slowly add 6mol/L NaOH solution\***, adjust the pH to **6.60 to 6.80\***, and stir for at least 15 minutes, at which time the solution should be gradually clear.  
*\*可按照比例折算使用其他浓度的氢氧化钠溶液, 但不建议使用氢氧化钠粉末, 因为局部过高的 pH 可能会产生损伤。建议缓慢调节溶液 pH 值, 避免超出此范围。*

*\*Other concentrations of NaOH solution can be converted in proportion. NaOH powder is not recommended in case local high pH causes damage. The solution pH is recommended to be adjusted slowly to avoid exceeding the pH range.*

5. **缓慢**加入 **2.20g/L** 碳酸氢钠粉末\*, 继续搅拌不少于 15 分钟。  
**Slowly add 2.20g/L sodium bicarbonate powder\*** and stir for at least 15 minutes.  
*\*碳酸氢钠粉末溶解后培养基 pH 会随着搅拌时间逐渐上升, 搅拌时间不宜过长。*  
*\*The pH will gradually increase by the mixing time after dissolution of sodium bicarbonate powder, and the mixing time should be controlled reasonably.*
6. 搅拌时间结束后, 检测培养基溶液 pH。若溶液 pH 低于 **7.00**, 则**缓慢**加入 6mol/L 氢氧化钠溶液调整 pH 至 **7.00~7.30** 范围\*。  
After stirring, test the pH of the medium solution. If the solution pH is below **7.00**, **slowly add 6mol/L NaOH solution** to adjust the pH to the range of **7.00 to 7.30\***.  
*\*pH 调整时应注意控制氢氧化钠溶液的使用量, 避免过度调整导致培养基最终渗透压超出标准范围。*  
*\*The addition amount of NaOH solution shall be controlled to adjust the pH, avoiding out-of-specification osmolality caused by excessive adjustment.*
7. 根据配制开始前确定的定容方式, 使用高纯度水定容至终体积或终重量\*, 继续搅拌 **10~15** 分钟。  
Add HPW to the final volume or final weight\* according to the determined q.s. method before preparation and stir for **10 to 15 minutes**.  
*\*培养基定容的具体操作和注意事项详见本说明书中“定容方式”部分的相关描述。*  
*\*Specific q.s. operations and precautions are described in the “Q.S. Method” section herein.*
8. 此时检测渗透压应为 **310~370mOsmol/kg**, 浊度应 **<3.0NTU**。  
At this time, expected osmolality is **310 to 370mOsmol/kg**, and expected turbidity is **<3.0NTU**.
9. 使用 **0.22µm 或 0.2µm** 除菌滤器\*过滤除菌。  
Sterilize with **0.22µm or 0.2µm sterile filter\***.  
*\*建议使用低结合类型滤膜, 如聚偏二氟乙烯 PVDF、聚醚砜 PES 或醋酸纤维素 CA。*  
*\*It is recommended to use a low-binding filter membrane type, such as PVDF, PES or CA.*
10. 立即使用, 或密封后在 **2 至 8°C** 条件下**避光保**

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Use the prepared medium immediately or seal and store at 2 to 8°C away from light.

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#### Precautions

1. 无需额外进行补料操作。  
Additional feed is not required.
2. DurePro® HD 灌流培养基含糖量 10.0g/L，灌流期间按需额外补充葡萄糖溶液。  
DurePro® HD Perfusion Medium contains 10.0g/L glucose, and additional glucose solution is supplemented as needed during perfusion.

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\*Except for process changes, other changes will not be notified.

©2026 CelluPro, Ltd. All rights reserved including data, graphics and images.  
**DurePro, CelluPro** and **Helix logo** are registered trademarks of CelluPro, Ltd.

---

烟台赛普生物技术有限公司  
CelluPro, Ltd.  
中国（山东）自由贸易试验区烟台片区烟台开发区北京中路 60 号  
60 Beijing Middle Road, Economic & Technological Development Area, Yantai,  
Shandong, China

电话 / Phone: +86 535 3573398  
电子邮箱 / E-mail: support@cellupro.cn  
www.cellupro.com